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GB

The Umbilical Hernia Support provides support in case of abdominal tissue infirmity in umbilical area. Supporting splints in dorsal area enhance the wearing comfort and avoid curling-up. The pad serves for re-positioning (pushing back) of the hernia. For this purpose the pad is reinforced by an inserted plastic core, and can be freely positioned on the entire support thanks to the hook strap.

Instructions for use:

Position pad onto the support by Velcro fixing and place onto the hernia port. Put support around the body and close it.

Indications:

Umbilical hernia, reducible umbilical hernias.

Contraindications:

Peripheral arterial disease (PAD), lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid, sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body.

Composition:

66 % Polyamide, 22 % polyester, 12 % spandex

Important instructions:

No reutilisation – this medical device is meant for the care of one patient only.
In case of doctor's order to wear support during night, avoid hindering blood circulation. In case of numbness loosen support and remove it if necessary. If complaints continue, please contact the doctor.

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